

Copper base, Type: Brass - Leaded, Disc 40 mm x 30 mm

Art. ID IMN-WG1
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,096	%		
	Bismuth (Bi)	[7440-69-9]	0,0013	%		
	Copper (Cu)	[7440-50-8]	60,99	%		
	Iron (Fe)	[7439-89-6]	0,0084	%		
	Manganese (Mn)	[7439-96-5]	0,16	%		
	Nickel (Ni)	[7440-02-0]	0,2	%		
	Phosphorus (P)	[7723-14-0]	0,029	%		
	Lead (Pb)	[7439-92-1]	0,71	%		
	Antimony (Sb)	[7440-36-0]	~0,062	%		
	Silicon (Si)	[7440-21-3]	~0,0046	%		
	Tin (Sn)	[7440-31-5]	0,29	%		
	Zinc (Zn)	[7440-66-6]	rem	%		