

Copper base, Type: Brass - Leaded, Disc 40 mm x 30 mm

Art. ID IMN-WG3
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,041	%		
	Bismuth (Bi)	[7440-69-9]	0,0057	%		
	Copper (Cu)	[7440-50-8]	58,2	%		
	Iron (Fe)	[7439-89-6]	0,31	%		
	Manganese (Mn)	[7439-96-5]	0,037	%		
	Nickel (Ni)	[7440-02-0]	0,029	%		
	Phosphorus (P)	[7723-14-0]	0,013	%		
	Lead (Pb)	[7439-92-1]	2,29	%		
	Antimony (Sb)	[7440-36-0]	0,018	%		
	Silicon (Si)	[7440-21-3]	~0,014	%		
	Tin (Sn)	[7440-31-5]	0,091	%		
	Zinc (Zn)	[7440-66-6]	rem	%		