

******WG4 - sold out*** Copper base, Type: Brass - Leaded, CRM, Set IMN-WG1 - IMN-WG6 only, see IMN-WG1-WG6, Disc 40 mm x 30 mm**

Art. ID IMN-WG4
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,073 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,014 | % | | |
| | Copper (Cu) | [7440-50-8] | 60,05 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,1 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,12 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,16 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,02 | % | | |
| | Lead (Pb) | [7439-92-1] | 1,41 | % | | |
| | Antimony (Sb) | [7440-36-0] | ~0,042 | % | | |
| | Silicon (Si) | [7440-21-3] | ~0,016 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,21 | % | | |
| | Zinc (Zn) | [7440-66-6] | rem | % | | |