

**\*\*\*\*WG4 - sold out\*\*\* Copper base, Type: Brass - Leaded, CRM, Set IMN-WG1 - IMN-WG6 only, see IMN-WG1-WG6, Disc 40 mm x 30 mm**

Art. ID                      IMN-WG4  
Unit                         disc  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,073	%		
	Bismuth (Bi)	[7440-69-9]	0,014	%		
	Copper (Cu)	[7440-50-8]	60,05	%		
	Iron (Fe)	[7439-89-6]	0,1	%		
	Manganese (Mn)	[7439-96-5]	0,12	%		
	Nickel (Ni)	[7440-02-0]	0,16	%		
	Phosphorus (P)	[7723-14-0]	0,02	%		
	Lead (Pb)	[7439-92-1]	1,41	%		
	Antimony (Sb)	[7440-36-0]	~0,042	%		
	Silicon (Si)	[7440-21-3]	~0,016	%		
	Tin (Sn)	[7440-31-5]	0,21	%		
	Zinc (Zn)	[7440-66-6]	rem	%		