

Copper base, Type: Brass - Leaded, Disc 40 mm x 30 mm

Art. ID IMN-WG5
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,058	%		
	Bismuth (Bi)	[7440-69-9]	0,0094	%		
	Copper (Cu)	[7440-50-8]	59,32	%		
	Iron (Fe)	[7439-89-6]	0,18	%		
	Manganese (Mn)	[7439-96-5]	0,074	%		
	Nickel (Ni)	[7440-02-0]	0,078	%		
	Phosphorus (P)	[7723-14-0]	0,016	%		
	Lead (Pb)	[7439-92-1]	1,66	%		
	Antimony (Sb)	[7440-36-0]	0,034	%		
	Silicon (Si)	[7440-21-3]	~0,022	%		
	Tin (Sn)	[7440-31-5]	0,14	%		
	Zinc (Zn)	[7440-66-6]	rem	%		