

Copper base, Type: Brass - Nickel, CRM, Disc 40 mm x 30 mm

Art. ID IMN-WP5
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,049 | % | | |
| | Arsenic (As) | [7440-38-2] | 0,021 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,021 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 0,026 | % | | |
| | Copper (Cu) | [7440-50-8] | 58,7 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,026 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,49 | % | | |
| | Nickel (Ni) | [7440-02-0] | 15,63 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,0027 | % | | |
| | Lead (Pb) | [7439-92-1] | ~1,8 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,028 | % | | |
| | Silicon (Si) | [7440-21-3] | ~0,03 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,33 | % | | |
| | Zinc (Zn) | [7440-66-6] | rem | % | | |