

Copper base, Type: Bronze, 85Cu5Sn5Pb5Zn, CRM, Chips 80 g

Art. ID IPT-10B
Unit 80 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Arsenic (As) | [7440-38-2] | 0,019 | % | | |
| | Copper (Cu) | [7440-50-8] | 85,2 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,211 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,33 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,003 | % | | |
| | Lead (Pb) | [7439-92-1] | 4,74 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,068 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,114 | % | | |
| | Tin (Sn) | [7440-31-5] | 4,61 | % | | |
| | Zinc (Zn) | [7440-66-6] | 4,73 | % | | |