

**Copper base, Type: Brass 2,45Pb, CRM, Chips 100 g**

Art. ID IPT-40  
Unit 100 g  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,002	%		
	Aluminium (Al)	[7429-90-5]	0,01	%		
	Cadmium (Cd)	[7440-43-9]	0,049	%		
	Copper (Cu)	[7440-50-8]	58,1	%		
	Iron (Fe)	[7439-89-6]	0,007	%		
	Nickel (Ni)	[7440-02-0]	0,0012	%		
	Lead (Pb)	[7439-92-1]	2,45	%		
	Antimony (Sb)	[7440-36-0]	0,023	%		
	Tin (Sn)	[7440-31-5]	0,18	%		
	Zinc (Zn)	[7440-66-6]	39,1	%		