

**Copper base, Type: Pure Copper, CRM, Chips 50 g**

Art. ID IPT-64  
Unit 50 g  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,001	%		
	Aluminium (Al)	[7429-90-5]	(<0,0006)	%		
	Arsenic (As)	[7440-38-2]	~0,0002	%		
	Bismuth (Bi)	[7440-69-9]	(<0,0001)	%		
	Copper (Cu)	[7440-50-8]	99,98	%		
	Iron (Fe)	[7439-89-6]	0,00045	%		
	Nickel (Ni)	[7440-02-0]	0,00018	%		
	Lead (Pb)	[7439-92-1]	0,00006	%		
	Antimony (Sb)	[7440-36-0]	~0,0002	%		
	Selenium (Se)	[7782-49-2]	(<0,0002)	%		
	Tin (Sn)	[7440-31-5]	(<0,0005)	%		
	Tellurium (Te)	[13494-80-9]	(<0,0001)	%		
	Zinc (Zn)	[7440-66-6]	~0,0001	%		