

Copper base, Type: Brass - Manganese, RM, Typical Analysis, Disc 45 mm x 25 mm

Art. ID PB-MS2
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,15	%		
	Bismuth (Bi)	[7440-69-9]	~0,031	%		
	Copper (Cu)	[7440-50-8]	64,44	%		
	Iron (Fe)	[7439-89-6]	0,063	%		
	Manganese (Mn)	[7439-96-5]	1	%		
	Nickel (Ni)	[7440-02-0]	0,34	%		
	Lead (Pb)	[7439-92-1]	0,64	%		
	Silicon (Si)	[7440-21-3]	0,079	%		
	Tin (Sn)	[7440-31-5]	0,64	%		