

**CHROMABOND QuEChERS Mix I citrate extraction mix consisting of: - 4 g MgSO<sub>4</sub> - 1 g NaCl - 0.5 g Na<sub>2</sub>H citrate x 1.5 H<sub>2</sub>O - 1 g Na<sub>3</sub> citrate x 2 H<sub>2</sub>O individually weighed and filled into 15 mL centrifuge tubes (PP) with screw cap (PE) pack of 50**

Art. ID	MN-730970
Unit	pack of 50
Deliverydetails	No Dangerous Good /not restricted