

VarioPrep HPLC guard column cartridge (short) VP 10/16 NUCLEODUR - polar universal guard column cartridge - length: 10 mm, ID: 16 mm pack of 2 - compatible with VarioPrep HPLC column with ID: 16 mm or ID: 21 mm - for all polar NUCLEODUR- and NUCLEOSIL phases - requires VP guard column holder 16 mm (short version, REF: 718256)

| | |
|-----------------|-----------------------------------|
| Art. ID | MN-762094.160 |
| Unit | PK |
| Deliverydetails | No Dangerous Good /not restricted |