

**CONOSTAN - Tungsten in 20 cSt blank oil, 5000 ppm**

Art. ID                      CON-150-500-745  
Unit                         50 g  
Deliverydetails            No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Tungsten (W)      | [7440-33-7] | 0,5                 | %    |        |        |
|                  | Tungsten (W)      | [7440-33-7] | 5000                | µg/g |        |        |