

Mixture of 26 components 300 µg/mL in Acetonitrile

Art. ID LM24-CUS-70123
Unit 6 x 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Aminobiphenyl	[92-67-1]	300	µg/mL		
	4,4'-Benzidine	[92-87-5]	300	µg/mL		
	4-Chloro-2-methylanilin e	[95-69-2]	300	µg/mL		
	2-Aminonaphthalene	[91-59-8]	300	µg/mL		
	4-Amino-2',3-dimethylaz obenzene	[97-56-3]	300	µg/mL		
	2-Amino-4-nitrotoluene (2-Methyl-5-nitroanilin e)	[99-55-8]	300	µg/mL		
	4-Chloroaniline	[106-47-8]	300	µg/mL		
	4-Methoxy-1,3-phenylene diamine	[615-05-4]	300	µg/mL		
	Bis-(4-aminophenyl)-met hane (4,4'-Methylenedia niline)	[101-77-9]	300	µg/mL		
	3,3'-Dichlorobenzidine	[91-94-1]	300	µg/mL		
	3,3'-Dimethoxybenzidine	[119-90-4]	300	µg/mL		
	3,3'-Dimethylbenzidine (o-Tolidine)	[119-93-7]	300	µg/mL		
	4,4'-Diamino-3,3'-dimet hyldiphenyl methane	[838-88-0]	300	µg/mL		
	2-Methoxy-5-methylanili ne	[120-71-8]	300	µg/mL		
	4,4'-Methylene-bis(2-ch loroaniline)	[101-14-4]	300	µg/mL		
	4-Aminophenylether (4,4 '-Oxydianiline)	[101-80-4]	300	µg/mL		
	4-Aminophenylthioether (4,4'-Thiodianiline)	[139-65-1]	300	µg/mL		
	2-Amino-1-methylbenzene	[95-53-4]	300	µg/mL		

(o-Toluidine)				
	2,4-Diaminotoluene	[95-80-7]	300	µg/mL
	2,4,5-Trimethylaniline	[137-17-7]	300	µg/mL
o-Anisidine	2-Anisidine	[90-04-0]	300	µg/mL
	2,4-Dimethylaniline	[95-68-1]	300	µg/mL
	2,6-Dimethylaniline	[87-62-7]	300	µg/mL
	4-Aminoazobenzene	[60-09-3]	300	µg/mL
	1,4-Phenylenediamine	[106-50-3]	300	µg/mL
	Aniline	[62-53-3]	300	µg/mL