

Reagecon GY1 to GY7 Colour Reference Set according to European Pharmacopoeia (EP)

Art. ID	REA-EPRSGYAMPSET
Unit	8 x 10 mL Set
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

Description

Reagecon's GY1 to GY7 Colour Reference Set is formulated as a set of reference solutions for the examination of the degree of colour of liquids in accordance with the European Pharmacopoeia Chapter 2.2.2. These products are produced by combining Reagecon Primary Colour Solutions.