

**Reagecon Y1 to Y7 Colour Reference Set according to European Pharmacopoeia (EP)**

Art. ID	REA-EPRSYAMPSET
Unit	8 x 10 mL Set
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

**Description**

Reagecon's Y1 to Y7 Colour Reference Set is formulated as a set of reference solutions for the examination of the degree of colour of liquids in accordance with the European Pharmacopoeia Chapter 2.2.2. These products are produced by combining Reagecon Primary Colour Solutions.