

Reagecon BY3 Colour Reference Solution according to European Pharmacopoeia (EP)

Art. ID	REA-EPBY103
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Reagecon's BY3 Colour Reference Set is formulated as a set of reference solutions for the examination of the degree of colour of liquids in accordance with the European Pharmacopoeia Chapter 2.2.2. These products are produced by combining Reagecon Primary Colour Solutions.