

## **Reagecon GY5 Colour Reference Solution according to European Pharmacopoeia (EP)**

Art. ID	REA-EPGY501
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

### Description

Reagecon's GY5 Colour Reference Solution is formulated as a reference solution for the examination of the degree of colouration of liquids in accordance with the European Pharmacopoeia Chapter 2.2.2. This product is produced by combining Reagecon Primary Colour Solutions.