

## **Tin - 1000 mg/l, 125 ml**

Art. ID ERA-037

Unit 125 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 III E1

### Description

One 125 mL standard containing Tin at 1000 mg/L in H<sub>2</sub>O with HCl. This CRM is covered by ERA's ISO Guide 34 accreditation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	1000	mg/L		