

Complex Cyanide - 1000 mg/l, 125 ml

Art. ID ERA-049
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1824 III E1

Description

One 125 mL standard containing Complex Cyanide at 1000 mg/L in H₂O with NaOH. This CRM is covered by ERA's ISO Guide 34 accreditation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Cyanide	Cyanide		1000	mg/L		