

Nutrients QC plus QC Standard

Art. ID	ERA-4023
Unit	each
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

Description

Two screw-cap vials yield up to 2 liters each after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Ammonia as N	Ammonia (as N)		0,25 - 10	mg/L		
Nitrate as N	Nitrate as N		0,25 - 10	mg/L		
ortho-Phosphate as P	Ortho-Phosphate as P		0,05 - 10	mg/L		
Phosphate as P	Phosphate as P		0,05 - 10	mg/L		
Total Kjeldahl Nitrogen	Total Kjeldahl Nitrogen (TKN)		0,25 - 10	mg/L		
Total phosphorus as P	Total Phosphorus as P		0,1 - 10	mg/L		