

Nutrients QC plus QC Standard

Art. ID ERA-4023
Unit each
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

Description

Two screw-cap vials yield up to 2 liters each after dilution.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|-------------------------|----------------------------------|------------|---------------------|------|--------|--------|
| Ammonia as N | Ammonia (as N) | | 0,25 - 10 | mg/L | | |
| Nitrate as N | Nitrate as N | | 0,25 - 10 | mg/L | | |
| ortho-Phosphate as P | Ortho-Phosphate as P | | 0,05 - 10 | mg/L | | |
| Phosphate as P | Phosphate as P | | 0,05 - 10 | mg/L | | |
| Total Kjeldahl Nitrogen | Total Kjeldahl Nitrogen (TKN) | | 0,25 - 10 | mg/L | | |
| Total phosphorus as P | Total Phosphorus as P | | 0,1 - 10 | mg/L | | |