

## **Solids QC plus QC Standard**

Art. ID	ERA-4033
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One screw-cap vial with a powder yields 1 liter after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Non-Filterable Residue (TSS)			mg/L			
Total Dissolved Solids at 180°C			mg/L			
Total Solids at 105°C			mg/L			
Total Suspended Solids			mg/L			