

Flame AA/ICP Cations

Art. ID	ERA-530
Unit	15 mL
Deliverydetails	Not yet checked

Description

One 15 mL screw-cap vial, unpreserved, yields up to 250 mL after dilution. Use with ICP, IC and AA methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	10-200	mg/L		
	Magnesium (Mg)	[7439-95-4]	10-200	mg/L		
	Potassium (K)	[7440-09-7]	5-100	mg/L		
	Sodium (Na)	[7440-23-5]	10-250	mg/L		