

Proficiency Testing - EDB/DBCP/TCP

| | |
|-----------------|-----------------|
| Art. ID | ERA-562 |
| Unit | 2 mL |
| Deliverydetails | Not yet checked |

Description

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Use with EPA method 8011. Each lot contains all analytes at .2-2.0 µg/L.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------------|------------|---------------------|------|--------|--------|
| | 1,2,3-Trichloropropane (TCP) | | 0,2-2 | µg/L | | |
| | 1,2-Dibromo-3-chloropropane (DBCP) | [96-12-8] | 0,2-2 | µg/L | | |
| | 1,2-Dibromoethane (EDB) | [106-93-4] | 0,2-2 | µg/L | | |