

Certified Reference Material - EDB/DBCP/TCP, WasteWatR?

Art. ID	ERA-692
Unit	2 mL
Deliverydetails	Not yet checked

Description

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Use with EPA method 8011. Each lot contains all analytes at .2 - 2.0 µg/L.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,3-Trichloropropane (TCP)		0,2-2	µg/L		
	1,2-Dibromo-3-chloropropane (DBCP)	[96-12-8]	0,2-2	µg/L		
	1,2-Dibromoethane (EDB)	[106-93-4]	0,2-2	µg/L		