

Ion Chromatography

Art. ID	ERA-981
Unit	15 mL
Deliverydetails	Not yet checked

Description

One 15 mL screw-cap vial yields up to 200 mL after dilution. Designed to calibrate or verify IC calibrations. Call for anion standards at lower levels.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bromide		0,2-20	mg/L		
	Chloride		0,2-20	mg/L		
	Fluoride		0,1-10	mg/L		
	Nitrate as N		0,2-20	mg/L		
	ortho-Phosphate as P		0,5-30	mg/L		
	Sulfate		0,5-30	mg/L		