

## **Free Cyanide - 1000 mg/l, 500 ml**

Art. ID ERA-997  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1824 III E1

### Description

One 500 mL standard containing Cyanide (free) at 1000 mg/L in H<sub>2</sub>O with NaOH. This CRM is covered by ERA's ISO Guide 34 accreditation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Free cyanide		1000	mg/L		