

Bismuth - 1000 mg/l, 125 ml

Art. ID ERA-K01
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

One 125 mL standard containing Bismuth at 1000 mg/L in H₂O with HNO₃. This CRM is covered by ERA's ISO Guide 34 accreditation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	1000	mg/L		