

Yttrium - 1000 mg/l, 125 ml

Art. ID ERA-K08
Unit 125 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

One 125 mL standard containing Yttrium at 1000 mg/L in H₂O with HNO₃. This CRM is covered by ERA's ISO Guide 34 accreditation.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1000	mg/L		
	Yttrium (Y)	[7440-65-5]	1000	mg/L		