

ICP-MS Major Cations

Art. ID ERA-TMS002
Unit 125 mL
Deliverydetails UN3264 III E1

Description

One 125 mL screw-cap poly bottle preserved with HNO₃.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Calcium (Ca) | [7440-70-2] | 50 | mg/L | | |
| | Magnesium (Mg) | [7439-95-4] | 50 | mg/L | | |
| | Potassium (K) | [7440-09-7] | 50 | mg/L | | |
| | Sodium (Na) | [7440-23-5] | 50 | mg/L | | |