

### **Custom Solution - 6 components in 2% nitric acid**

Art. ID LM24-CUS-70125  
Unit 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	50	mg/L		
	Lead (Pb)	[7439-92-1]	50	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	10	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	50	mg/L		