

### **Calibration Standard #1 Method 200,7, Varied conc, in 2-5% Nitric Acid**

Art. ID ACCU-M-200.7-01-1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	50	µg/mL		
	Arsenic (As)	[7440-38-2]	100	µg/mL		
	Barium (Ba)	[7440-39-3]	10	µg/mL		
	Boron (B)	[7440-42-8]	20	µg/mL		
	Cadmium (Cd)	[7440-43-9]	20	µg/mL		
	Calcium (Ca)	[7440-70-2]	100	µg/mL		
	Copper (Cu)	[7440-50-8]	20	µg/mL		
	Manganese (Mn)	[7439-96-5]	20	µg/mL		
	Selenium (Se)	[7782-49-2]	50	µg/mL		
	Silver (Ag)	[7440-22-4]	5	µg/mL		