

**Method 200,7 Calibration Standard 2, Varied conc, in 2-5% Nitric Acid, tr Hydrofluoric acid**

Art. ID ACCU-M-200.7-02R-1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium (Li)	[7439-93-2]	50	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	100	µg/mL		
	Potassium (K)	[7440-09-7]	200	µg/mL		
	Sodium (Na)	[7440-23-5]	100	µg/mL		
	Strontium (Sr)	[7440-24-6]	10	µg/mL		
	Titanium (Ti)	[7440-32-6]	100	µg/mL		