

Method 200,7 Calibration Standard 3, Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-M-200.7-03R-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cerium (Ce)	[7440-45-1]	20	µg/mL		
	Cobalt (Co)	[7440-48-4]	20	µg/mL		
	Phosphorus (P)	[7723-14-0]	100	µg/mL		
	Vanadium (V)	[7440-62-2]	20	µg/mL		