

**Calibration Standard #4 Method 200,7, Varied conc, in 2-5% Nitric Acid, tr Hydrofluoric acid**

Art. ID ACCU-M-200.7-04-5  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	µg/mL		
	Chromium (Cr)	[7440-47-3]	50	µg/mL		
	Silicon (Si)	[7440-21-3]	100	µg/mL		
	Tin (Sn)	[7440-31-5]	40	µg/mL		
	Zinc (Zn)	[7440-66-6]	50	µg/mL		