

**Calibration Standard #5 Method 200,7, Varied conc, in 2-5% Nitric Acid**

Art. ID ACCU-M-200.7-05-1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	10	µg/mL		
	Iron (Fe)	[7439-89-6]	100	µg/mL		
	Lead (Pb)	[7439-92-1]	100	µg/mL		
	Magnesium (Mg)	[7439-95-4]	100	µg/mL		
	Nickel (Ni)	[7440-02-0]	20	µg/mL		
	Thallium (Tl)	[7440-28-0]	50	µg/mL		