

Spiking Std #1 Method 200,7 (1991), Varied conc, in Water, tr Hydrofluoric acid

Art. ID ACCU-M-200.7-SP-01-R
Unit 50 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	400	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	200	µg/mL		
	Silicon (Si)	[7440-21-3]	2000	µg/mL		
	Phosphorus (P)	[7723-14-0]	400	µg/mL		