

Spiking Std #5 Method 200,7 (1991), Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-M-200.7-SP-05-R
Unit 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	200	µg/mL		
	Cadmium (Cd)	[7440-43-9]	100	µg/mL		
	Lead (Pb)	[7439-92-1]	200	µg/mL		
	Selenium (Se)	[7782-49-2]	400	µg/mL		
	Thallium (Tl)	[7440-28-0]	400	µg/mL		