

### **Azspiracid-3 standard solution in Methanol**

Art. ID	NRC-CRM-AZA3-b
Unit	0,5 mL
Deliverydetails	De Minimis UN1230 (6.1) II E2

#### Description

Each ampoule contains approximately 0.5 mL of solution with  $1.25 \pm 0.05 \mu\text{mol/L}$  AZA3 in methanol. It is designed for use as an instrument calibration solution to aid the analyst in the determination of AZA3 in plankton and shellfish samples.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Azspiracid-3	[265996-93-8]	1,43	$\mu\text{mol/L}$		