

Trace elements in water (lot 0622)

Art. ID NW-TMDA-61.3
Unit 250 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	58,6 ± 6,4	µg/L		
	Arsenic (As)	[7440-38-2]	35,2 ± 3,5	µg/L		
	Barium (Ba)	[7440-39-3]	64,3 ± 4	µg/L		
	Beryllium (Be)	[7440-41-7]	37 ± 3,5	µg/L		
	Bismuth (Bi)	[7440-69-9]	19,9 ± 4,2	µg/L		
	Boron (B)	[7440-42-8]	48,7 ± 6,6	µg/L		
	Cadmium (Cd)	[7440-43-9]	59,3 ± 4	µg/L		
	Chromium (Cr)	[7440-47-3]	67,7 ± 4,9	µg/L		
	Cobalt (Co)	[7440-48-4]	63,7 ± 4,4	µg/L		
	Copper (Cu)	[7440-50-8]	62,1 ± 4,9	µg/L		
	Gallium (Ga)	[7440-55-3]	8,42 ± 0,5	µg/L		
	Iron (Fe)	[7439-89-6]	79,7 ± 6,6	µg/L		
	Lead (Pb)	[7439-92-1]	56,1 ± 4,8	µg/L		
	Lithium (Li)	[7439-93-2]	34,8 ± 3,5	µg/L		
	Manganese (Mn)	[7439-96-5]	76,4 ± 5,1	µg/L		
	Molybdenum (Mo)	[7439-98-7]	73,1 ± 5,4	µg/L		
	Nickel (Ni)	[7440-02-0]	56,9 ± 4,3	µg/L		
	Rubidium (Rb)	[7440-17-7]	3,23 ± 0,24	µg/L		
	Selenium (Se)	[7782-49-2]	40,2 ± 4,3	µg/L		
	Strontium (Sr)	[7440-24-6]	114 ± 7	µg/L		
	Thallium (Tl)	[7440-28-0]	37,3 ± 3,2	µg/L		
	Tin (Sn)	[7440-31-5]	56 ± 5,4	µg/L		
	Titanium (Ti)	[7440-32-6]	37,3 ± 3	µg/L		
	Tungsten (W)	[7440-33-7]	2,07 ± 0,21	µg/L		
	Uranium (U)	[7440-61-1]	36,3 ± 2,5	µg/L		
	Vanadium (V)	[7440-62-2]	70,5 ± 5	µg/L		
	Zinc (Zn)	[7440-66-6]	72,8 ± 6,6	µg/L		