

### **Custom Solution - 6 components in 2% nitric acid**

Art. ID LM24-CUS-70133  
Unit 30 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1760 III E1

#### Description

Stored under normal laboratory conditions. Conservation 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium (Li)	[7439-93-2]	10	µg/mL		
	Magnesium (Mg)	[7439-95-4]	10	µg/mL		
	Yttrium (Y)	[7440-65-5]	10	µg/mL		
	Cerium (Ce)	[7440-45-1]	10	µg/mL		
	Thallium (Tl)	[7440-28-0]	10	µg/mL		
	Cobalt (Co)	[7440-48-4]	10	µg/mL		