

## **PremiSolv™ ICP solvent**

Art. ID	CON-150-700-002
Unit	18,9 L
Deliverydetails	No Dangerous Good /not restricted

### Description

PremiSolv ICP Solvent is a replacement for Low Odor Kerosene and Xylene for use as a diluent in ICP/DCP analysis of metals in oil and other organic fluids. PremiSolv ICP solvent is supplied in one gallon containers for shipment via common carrier or drum size quantities for high volume users. PremiSolv ICP solvent features: 1. Supplied with Certificate of Analysis listing the concentration of 34 different metals \_\_ 2. Extremely low toxicity compared to Kerosene or Xylene \_\_ 3. No Sulfur/Aromatic odor for increased worker comfort and less exposure \_\_ 4. Non-Hazardous for shipping/disposal