

PORCINE pDNA CALIBRANT

Art. ID ERM-AD483
Unit set of 3 vials
Deliverydetails Dry ice shipment

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ERM-AD483a: Copy number concentration of the p lasmid			126 ± 18	cp/μL		
ERM-AD483b: Copy number concentration of the p lasmid			34 ± 6	cp/μL		
ERM-AD483c: Copy number concentration of the p lasmid			9 ± 3	cp/μL		