

Cu (1.0 mm wire)

Art. ID	IRMM-522D
Unit	1 m (7.0 g)
Deliverydetails	No Dangerous Good /not restricted

Description

IRMM-522D consists of copper with a nominal purity of $\geq 99.995\%$. The isotopic abundances correspond to the natural composition. Cobalt and silver, present as impurities and causing interference during application, are certified. The material is available as 0.5 mm and 1 mm diameter wire (1 m units), 1 mm thick foil (20 cm² units) and 0.1 mm thick foil (100 cm² units).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cobalt (Co)	[7440-48-4]	<0,05	mg/kg		
	Silver (Ag)	[7440-22-4]	0,95 ± 0,04	mg/kg		