

Custom Solutions - 9 components in 7% HNO₃

Art. ID LM24-CUS-70137

Unit 2x 500 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	40	mg/L		
	Manganese (Mn)	[7439-96-5]	25	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	1000	mg/L		
	Magnesium (Mg)	[7439-95-4]	500	mg/L		
	Sodium (Na)	[7440-23-5]	500	mg/L		
	Potassium (K)	[7440-09-7]	200	mg/L		
	Strontium (Sr)	[7440-24-6]	10	mg/L		
	Lithium (Li)	[7439-93-2]	10	mg/L		