

Custom Solution - 3 components in 2% HNO₃

Art. ID LM24-CUS-70138
Unit 2x 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	200	mg/L		
	Manganese (Mn)	[7439-96-5]	25	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		