

Diesel - CFPP - Cold Filter Plugging Point (Main Method: ASTM D 6371 - EN 116)

Art. ID ROF-CF-GO-333
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------------|------------|---------------------|------|--------|--------|
| | Cold filter plugging point (CFPP) | | -16 | °C | | |