

**Test Dust (0,09 - 0,20 mm) (according to DIN EN 60312:2008-11)**

Art. ID	NOS-006-7.5KG
Unit	can (7,5 kg)
Deliverydetails	No Dangerous Good /not restricted

**Description**

Test dust is a product to evaluate the performance characteristics of domestic used vacuum cleaner. The sand is composed according to DIN EN 60312:2008-11 chapter 5.1.2.2. It is sand with a clearly defined grain size distribution in the range of 0.09 mm – 0.20 mm. To determine the performance of a vacuum cleaner the sand is spread on and rolled into a standardized carpet. The performance of the vacuum cleaner is the ratio between the hoovered mass and the total mass.