

15(E)-Tetracosenoic acid, C24:1, Purity: >98%

Art. ID LA-10-2408-4
Unit 5 mg

Description

Synonyms: trans-Nervonic acid / trans-15-Tetracosenoic acid / 15-trans-Tetracosenoic acid

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------|---------------|---------------------|------|--------|--------|
| | 15(E)-Tetracosenoic aci d | [115863-91-7] | | | | |